
Grip-Lok®

REF 3300MW-TA-3

- EN Universal Fixation
- DE Universalbefixierung
- FR Fixation universelle
- IT Dispositivo di fissaggio universale
- NL Universele fixatie
- SV Universell fixering
- ES Fijación universal
- PT Fixação Universal
- EL Καθήλωση γενικής χρήσης
- FI Yleiskäyttöinen kiinnitys
- DA Universalfixering
- IS Alhliða festing
- NO Universell fiksering
- PL Mocowanie uniwersalne
- CS Univerzální fixace
- SK Univerzálna fixácia

44U00023 REV 01

TeleflexMedical Europe
IDA Business and Technology Park
Dublin Road, Athlone, Co. Westmeath, Ireland
Phone + 353 (0)9 06 46 08 00
www.teleflexmedical.com

Teleflex®

TIDI Products, LLC
570 Enterprise Drive
Neenah, WI 54956, USA
Phone: + 1.920.751.4300
tidiproducts.com/securement

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UK CA 0086
UK Responsible Person:
Emergo Consulting (UK) Limited
c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

Grip-Lok is a registered trademark of TIDI Products, LLC.
May be patented: see www.tidiproducts.com/virtual-patent-marking

2023-01-03

INSTRUCTIONS

EN

Intended for use in the general population. English

- Select the area for the placement of the Grip-Lok. Note: catheter, tube or line should be already placed/inserted.
- Prepare the skin according to the standard hospital protocol for dressing application. Skin prep or hair removal may be required on some patients for better adhesion.
- Open the top flap and slide the Grip-Lok under the catheter, tube or line, centering it above the adhesive area.
- Pull the paper backing from one side of the Grip-Lok, then the other, to secure in desired position on the skin.
- Remove the interior liner to expose the adhesive areas.
- Place the catheter, tube or line in the center of the exposed adhesive.
- Secure top flap over the catheter, tube or line and apply gentle pressure to the center fabric area over the adhesive.

Direction: For additional skin protection, use skin prep pad prior to application of stabilization device.

Note: Use of an alcohol swab may help removal of securement device from the skin.

⚠ If you experience severe redness, itching, swelling, or irritation of the skin, consult your physician as this may be a sign of an allergic reaction.

⚠ Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

⚠ Replace securement device if soiled or saturated in fluid or if device shows signs of wear or damage.

Note: If being used in a medical facility, dispose according to medical facility and municipal guidelines.

⊘ Do not use if the product's sterile barrier system or its packaging is compromised.

Report any serious incident related to device use to TIDI Products and the Member State competent authority.

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