



# Safety **+**Plus GOWNS™

## P2® SafetyPlus™ GOWNS IN TWO NEW OPTIONS FOR COMFORT AND PROTECTION

Exposure to infectious diseases is at an all-time high and the use of personal protection equipment is critical to assist in reducing the spread of hospital-contracted infections.

### P2® SafetyPlus™ SMS GOWN

P2® SafetyPlus™ SMS Gowns are made from materials that are breathable and comfortable while still providing protection for the healthcare worker. Recommended for use when the potential of light to moderate fluids and/or contact isolation is anticipated.

- Yellow three-layer fabric
- Elastic cuff to prevent sleeve from moving up
- Closed back for full protection
- Ties at neck and waist for secure fit

10 bags of 10 gowns | 100 gowns/case

SIZE	PRODUCT NUMBER
------	----------------

X-Large	8579
---------	------

### P2® SafetyPlus™ POLY COATED GOWN

These gowns are made from a soft breathable material with an impervious coating to protect the healthcare worker. Recommended for use when the potential of moderate fluids and/or contact isolation is anticipated.

- White Polyethylene coated Polypropylene
- Stirrup cuff for wrist protection
- Open back for increased airflow
- Closed loop neck for easy donning and tear off
- Dispenser box is compact and uses 75% less space than conventional isolation gown packaging

5 dispensers of 8 gowns | 40 gowns/case

SIZE	PRODUCT NUMBER
------	----------------

X-Large	8578
---------	------



#8578

#8579

Personal protective equipment will be considered 'appropriate' only if it does not permit blood or other potentially infectious materials to pass through to or reach the employee's work clothes, street clothes, undergarments, skin, eyes, mouth, or other mucous membranes, under normal conditions of use and for the duration of time which the protective equipment will be used

– OSHA Standard 29  
CFR Part 1910.1030(d)(2)(i)

For full information on TIDI's P2® Impervious gown, refer to our P2® brochure or visit our website.

TIDI Products, LLC | 570 Enterprise Drive | Neenah, WI 54956  
Customer Service 800.521.1314 | FAX 800.837.7770 | [www.tidiproducs.com](http://www.tidiproducs.com)  
08-GWNFLIER-00

