

Posey® Fall Management Bracelets


Application Instructions



DESCRIPTION OF PRODUCT: Posey Fall Management Bracelets offer a non-intrusive way to clearly identify a high fall risk patient. Strong and durable, yet lightweight for patient comfort. Available in five colors. “Fall Risk” Bracelet available in yellow only.

Application Instructions

1. Size bracelet around patient’s wrist (Fig. 1).
2. Line up desired hole and snap fastener closed (Fig 2).
3. You should be able to easily insert one finger between the bracelet and the patient’s limb (Fig. 3).
4. Trim off any excess wristband material with scissors, leaving smooth, rounded corners.



Fall Management Bracelets

- | | |
|-------------------------|-----------------------------------|
| REF 6247R Red | REF 6247G Green |
| REF 6247O Orange | REF 6247Y Yellow |
| REF 6247P Purple | REF 6251Y Yellow Fall Risk |

TIDI Products, LLC • 570 Enterprise Drive, Neenah, WI 54956 USA
 Phone: 1.800.521.1314 • International: +1.920.751.4036 • www.tidiproducts.com

© 2024 TIDI Products, LLC. All rights reserved. Posey is a registered trademark of TIDI Products, LLC.

M6159 REV D 04-02-2024

Posey® Fall Management Bracelets


Application Instructions



DESCRIPTION OF PRODUCT: Posey Fall Management Bracelets offer a non-intrusive way to clearly identify a high fall risk patient. Strong and durable, yet lightweight for patient comfort. Available in five colors. “Fall Risk” Bracelet available in yellow only.

Application Instructions

1. Size bracelet around patient’s wrist (Fig. 1).
2. Line up desired hole and snap fastener closed (Fig 2).
3. You should be able to easily insert one finger between the bracelet and the patient’s limb (Fig. 3).
4. Trim off any excess wristband material with scissors, leaving smooth, rounded corners.



Fall Management Bracelets

- | | |
|-------------------------|-----------------------------------|
| REF 6247R Red | REF 6247G Green |
| REF 6247O Orange | REF 6247Y Yellow |
| REF 6247P Purple | REF 6251Y Yellow Fall Risk |

TIDI Products, LLC • 570 Enterprise Drive, Neenah, WI 54956 USA
 Phone: 1.800.521.1314 • International: +1.920.751.4036 • www.tidiproducts.com

© 2024 TIDI Products, LLC. All rights reserved. Posey is a registered trademark of TIDI Products, LLC.

M6159 REV D 04-02-2024