

Fall Management Blankets

Posey Fall Management Blankets provide a discreet way to identify a high fall risk patient. The color helps staff members quickly identify a high fall risk patient. Made of polyester. Meets CA #117, Section E Flame Retardancy Standards. Available in five colors; red, orange, purple, green and yellow.

- | | | | |
|----------------------------------|-------------------------|----------------------------------|-------------------------|
| REF 6248R Red (1 each) | 6248BR (36 each) | REF 6248O Orange (1 each) | 6248BO (36 each) |
| REF 6248P Purple (1 each) | 6248BP (36 each) | REF 6248G Green (1 each) | 6248BG (36 each) |
| REF 6248Y Yellow (1 each) | 6248BY (36 each) | | |

Application Instructions:

Posey Fall Management Blankets can be draped over a patient’s lap to preserve dignity while sitting in a chair, (Fig. 2) or simply hung over the foot of the bed while the patient sleeps (Fig. 1).

Washing Instructions:

Wash with similar colors.
Machine wash warm, tumble dry low.



Storage And Handling

This device is designed for use in normal indoor environments. This device may be stored in ambient warehouse temperatures at normal humidity levels. Avoid excess moisture or high humidity that may damage product materials.

Disposal

⚠️WARNING Properly dispose of the product per facility’s policy for BIOHAZARDOUS materials.



TIDI Products, LLC • 570 Enterprise Drive, Neenah, WI 54956 USA
Phone: 1.800.521.1314 • International: +1.920.751.4036 • www.tidiproducts.com

© 2024 TIDI Products, LLC. All rights reserved. Posey is a registered trademark of TIDI Products, LLC.

M6175 REV D 03-20-2024

Fall Management Blankets

Posey Fall Management Blankets provide a discreet way to identify a high fall risk patient. The color helps staff members quickly identify a high fall risk patient. Made of polyester. Meets CA #117, Section E Flame Retardancy Standards. Available in five colors; red, orange, purple, green and yellow.

- | | | | |
|----------------------------------|-------------------------|----------------------------------|-------------------------|
| REF 6248R Red (1 each) | 6248BR (36 each) | REF 6248O Orange (1 each) | 6248BO (36 each) |
| REF 6248P Purple (1 each) | 6248BP (36 each) | REF 6248G Green (1 each) | 6248BG (36 each) |
| REF 6248Y Yellow (1 each) | 6248BY (36 each) | | |

Application Instructions:

Posey Fall Management Blankets can be draped over a patient’s lap to preserve dignity while sitting in a chair, (Fig. 2) or simply hung over the foot of the bed while the patient sleeps (Fig. 1).

Washing Instructions:

Wash with similar colors.
Machine wash warm, tumble dry low.



Storage And Handling

This device is designed for use in normal indoor environments. This device may be stored in ambient warehouse temperatures at normal humidity levels. Avoid excess moisture or high humidity that may damage product materials.

Disposal

⚠️WARNING Properly dispose of the product per facility’s policy for BIOHAZARDOUS materials.



TIDI Products, LLC • 570 Enterprise Drive, Neenah, WI 54956 USA
Phone: 1.800.521.1314 • International: +1.920.751.4036 • www.tidiproducts.com

© 2024 TIDI Products, LLC. All rights reserved. Posey is a registered trademark of TIDI Products, LLC.

M6175 REV D 03-20-2024