

### Removable "ID" Bracelet


- Temporarily replaces permanent "ID" when permanent "ID" interferes with IV site.
- Allows "ID" to remain on patient at all times - not taped to outside of isolette.
- Easily moved from arm to leg without disrupting skin integrity.
- No adhesive "tape burns" when removed.
- Soft foam with hook and loop closure. No tape!
- One dozen per box.

**REF 4648** Newborn, 5½"L x 1¼"W (14 cm x 3 cm)


**REF 4649** Preemie, 2¾"L x 1¼"W (7 cm x 3 cm)

### Storage and Handling

This device is designed for use in normal indoor environments. This device may be stored in ambient warehouse temperatures at normal humidity levels. Avoid excess moisture or high humidity that may damage product materials.



4648



4649

### Disposal

**⚠ WARNING** Properly dispose of the product per facility's policy for BIOHAZARDOUS materials.

### Removable "ID" Bracelet


- Temporarily replaces permanent "ID" when permanent "ID" interferes with IV site.
- Allows "ID" to remain on patient at all times - not taped to outside of isolette.
- Easily moved from arm to leg without disrupting skin integrity.
- No adhesive "tape burns" when removed.
- Soft foam with hook and loop closure. No tape!
- One dozen per box.

**REF 4648** Newborn, 5½"L x 1¼"W (14 cm x 3 cm)


**REF 4649** Preemie, 2¾"L x 1¼"W (7 cm x 3 cm)

### Storage and Handling

This device is designed for use in normal indoor environments. This device may be stored in ambient warehouse temperatures at normal humidity levels. Avoid excess moisture or high humidity that may damage product materials.



4648



4649

### Disposal

**⚠ WARNING** Properly dispose of the product per facility's policy for BIOHAZARDOUS materials.